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1. Herman A. et al., Allergic contact dermatitis caused by isobornyl acrylate in Freestyle Libre, a new introduced glucose sensor, Contact Dermatitis, 2017;77(6):367-373

2. Korsgaard Berg A. et al., Skin Problems Associated with Insulin Pumps and Sensors in Adults with Type 1 Diabetes: A Cross Sectional Study, Diabetes Technol Ther, 2018;20(7):475-482

3. Svedman C. et al., Changes in adhesive ingredients in continuous glucose monitoring systems may induce new contact allergy pattern, Contact Dermatitis, 2021;84(6):439-446

4. Heinemann L, Kamann S. Adhesives Used for Diabetes Medical Devices: A Neglected Risk With Serious Consequences? J Diabetes Sci Technol. 2016;10(6):1211-1215. Published 2016 Aug 25. 9

5. Messer LH, Berget C, Beston C, Polsky S, Forlenza GP. Preserving Skin Integrity with Chronic Device Use in Diabetes. Diabetes Technol Ther. 2018;20(S2):S254-S264.

6. Lind M, Polonsky, W, Hirsch, I et al. Continuous Glucose Monitoring vs Conventional Therapy for Glycemic Control in Adults With Type 1 Diabetes Treated With Multiple Daily Injections – The GOLD Randomized Clinical Trial. JAMA. 2017;317(4):379-387.

7. Engler R, Routh TL, Lucisano JY. Adoption Barriers for Continuous Glucose Monitoring and Their Potential Reduction With a Fully Implanted System: Results From Patient Preference Surveys. Clin Diabetes. 2018;36(1):50-58

8. McNichol L, Lund C, Rosen T, Gray M., "Medical Adhesives and Patient Safety: State of the Science", J Wound Ostomy Continence Nurs. 2013;40(4):1-15

9. Christiansen MP, Kla LJ, Brazg R, et al., A Prospective Multicenter Evaluation of the Accuracy of a Novel Implanted Continuous Glucose Sensor. PRECISE II. Diabetes Technology & Therapeutics. 2018;20(3):197-206

10. Garg, S. K. et al., Evaluation of Accuracy and Safety of the Next-Generation Up to 180-Day Long-Term Implantable Eversense Continuous Glucose Monitoring System: The PROMISE Study. Diabetes Technology & Therapeutics, 2021;24(2), 1-9. DOI: 10.1089/dia.2021.0182

11. Barnard Katharine D, Kropff Jort, et al. Acceptability of Implantable Continuous Glucose Monitoring Sensor. Journal of Diabetes Science and Technology 2018;12(3):634-638.

Distribué par :
Ascensia Diabetes Care Switzerland AG
Peter Merian-Strasse 90
4052 Basel | Switzerland
061 544 79 90 | ascensia-diabetes.ch

Fabriqué par : Senseonics, Inc.
20451 Seneca Meadows Parkway
Germantown, MD 20876-7005 USA

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